

Shake-Rx

Automated Vial Shaking Device



Let it shake for you!

Shake-Rx is designed to accelerate the time-consuming vial shaking process after the drug reconstitution. Adjustable shaking speed and timer functions help dissolve up to 27 vials simultaneously.



Smart Solution for Vial Shaking Process



- Small footprint
- Plug & Play
- Multiple vial shaking (up to 27)
- Timer adjustment
- 30 speed levels
- 3 replaceable vial housing area
- Compatible with different vial types

Technical Information

Load Capacity:	For small vials: up to 27, for medium vials: up to 15, for big vials: up to 9
Timer range:	Starting from 15sec to 60min
Weight:	4.7 kg
Dimensions (w x d x h):	23x28x14 cm
Power:	12VDC Power Adaptor - 16W max.



KAPSAM HEALTH PRODUCTS LTD.

■ Cumhuriyet Mah. İsmet İnönü Cad. No: 57

41420 Çayırova, Kocaeli - Türkiye

☎ +90(262) 2908801 🌐 www.kapsamsaglik.com.tr

